



eurogine
Technology for Women's Healthcare

GOLD® T

Références
4267365 / 4267332 / 4267357

DESCRIPTION

Intra-uterine device, with a "Y" shaped structure. Can be placed in the insertion tube by pulling the wires.

Copper surface of 375 mm², with a gold core.

Intra-uterine contraception. Sterile. Product for single use.

EC mark, non-active implantable device in accordance with Directive 93/42/EEC Class III.



More for Less...



TECHNICAL FUNCTIONS

COMPONENT	COMPOSITION	DIMENSIONS			
IUD structure	Polyethylen + Barium Sulfate	Ref.	Size	Width	Length
		4267365	Normal	31 mm	33 mm
		4267332	Mini	24 mm	30.50 mm
		4267357	Maxi	36.50 mm	38 mm
Wire	Copper purity 99,99%	Diam.	Length	Area	Diam. Gold Core
		0.40 mm	302 mm	375 mm ²	0.10 mm
Withdrawal wire	Polyethylen	Diam.	Length		
		0.25 mm	380 mm x 2		
Insertion tube	Polyethylen	External Diam.			
		3.80 mm			
Insertion device	ABS	Hysterometry from 5 to 10 cm with 0.5 cm interval			

STERILISATION METHOD

Ethylen Oxide

INSERTION SYSTEM

Adjustable insertion system that ensures, by means of the "locking plunger", the exact placement of the IUD in the bottom of the uterus, reducing the risk of uterine puncture.

PRESENTATION

1 x 1 packaging

MAXIMUM DURATION "IN SITU"

5 years